



## CHEMOTHERAPY PROTECTIVE APRON WITH SLEEVES



**BioClean-C™** chemotherapy protective aprons with sleeves have been carefully researched, designed and constructed to provide the wearer with maximum protection against a range of chemotherapy drugs, having been tested using the new ASTM D 6978-05 standard, without compromising the comfort of the garment.

These chemotherapy protective aprons with sleeves have been manufactured from the revolutionary new **Clean-Tough™** fabric which is lightweight, low linting, hard-wearing and compatible for use in sterile and non-sterile cleanrooms and hospital pharmacies.

The unique design of our chemotherapy protective aprons with sleeves ensures maximum protection for the wearer. They are available in three sizes to suit all wearers and in sterile and non-sterile versions.



GB04/61092

NITRITEX LTD.	T: +44 (0) 1638 663338	E: <a href="mailto:info@nitritex.com">info@nitritex.com</a>
NITRITEX ASIA PTE LTD.	T: +65 6536 9123	E: <a href="mailto:nitritex@pacific.net.sg">nitritex@pacific.net.sg</a>
NITRITEX PUERTO RICO Inc.	T: +1 787 309 3979	E: <a href="mailto:info@nitritex.com.pr">info@nitritex.com.pr</a>
NITRITEX CANADA LTD.	T: +1 (905) 946 9539	E: <a href="mailto:info@nitritexcanada.com">info@nitritexcanada.com</a>
NITRITEX MALAYSIA Sdn. Bhd.	T: +(603) 5569 3857	E: <a href="mailto:info@nitritex.com.my">info@nitritex.com.my</a>



## CHEMOTHERAPY PROTECTIVE APRON WITH SLEEVES

Colour:	Blue.												
Features:	Adjustable neck fastening. Tie-tapes at waist. Attached sleeves with elasticated cuffs.												
Material:	Clean-Tough™ fabric.												
Printing:	Our special print process ensures that bags can be wiped down with cleaning or sterilising fluids without fear of cross contamination from printing inks.												
Sterilisation Method:	Gamma irradiation where applicable.												
Quality Systems:	Manufactured in a facility holding the registration: ISO 9001:2001.												
Sampling & Inspection:	In accordance with ISO 2859.												
Storage Instructions:	Store in a dry, cool place (<40°C) away from direct sunlight and fluorescent light.												
Bar Codes:	A unique bar code on the outer carton (ITF-14). This allows automatic data entry into Stock Control Systems.												
Lot Number:	This appears on each individual bag and the outer carton to provide traceability.												
Packaging:	One piece per sealed inner bag; One sealed inner bag per sealed outer bag; 40 outer bags per carton liner; One carton liner per carton (40 pieces).												
Compatibility:	ISO Class 4. EU GMP Grade A (Sterile Version).												
Re-order Codes and Sizes:	<table><tr><td>Non-sterile</td><td>Sterile</td><td>Size</td></tr><tr><td>BCAS-S</td><td>S-BCAS-S</td><td>Small</td></tr><tr><td>BCAS-M</td><td>S-BCAS-M</td><td>Medium</td></tr><tr><td>BCAS-L</td><td>S-BCAS-L</td><td>Large</td></tr></table>	Non-sterile	Sterile	Size	BCAS-S	S-BCAS-S	Small	BCAS-M	S-BCAS-M	Medium	BCAS-L	S-BCAS-L	Large
Non-sterile	Sterile	Size											
BCAS-S	S-BCAS-S	Small											
BCAS-M	S-BCAS-M	Medium											
BCAS-L	S-BCAS-L	Large											

