

Omega Neotech Sterile™

STERILE NON-LATEX EXAMINATION GLOVES



ONBS
Master Code

Omega Neotech™ sterile powderfree ambidextrous non-latex medical gloves mark an important advance in non-latex thin glove technology

These non-latex medical gloves are 100% latex-free, eliminating allergy problems normally associated with natural latex proteins.

The gloves have excellent chemical resistance and the soft formulation provides superior comfort and feel. They are available in a wide range of sizes to suit all users.



GB04/61092




NITRITEX LTD.	T: +44 (0) 1638 663338	E: info@nitritex.com
NITRITEX ASIA PTE LTD.	T: +65 6536 9123	E: nitritex@pacific.net.sg
NITRITEX PUERTO RICO Inc.	T: +1 787 309 3979	E: info@nitritex.com.pr
NITRITEX CANADA LTD.	T: +1 (905) 946 9539	E: info@nitritexcanada.com
NITRITEX MALAYSIA Sdn. Bhd.	T: +(603) 5569 3857	E: info@nitritex.com.my

 **Nitritex**
www.nitritex.com

Omega Neotech Sterile™

STERILE NON-LATEX EXAMINATION GLOVES

Colour:	Blue.	 Master Code
Length:	240mm (9.5").	
Material:	Polychloroprene.	
Type:	Ambidextrous.	
Surface:	Finger textured.	
Other Features:	Beaded cuff, powderfree, latex-free.	
Protein Level:	Polychloroprene contains no natural latex proteins.	
Sterilisation Method:	Gamma Irradiation.	
Packing:	Boxed.	
AQL:	1.5 for holes (1000ml water test).	
Physical Properties:	Manufactured to exceed the requirements of the European Medical Gloves Standard EN455 (Parts 1, 2 & 3).	
Quality Systems:	Manufactured in a facility holding the ISO 9001:2000 registration and Directive 93/42/EEC Annex V.	
Sampling & Inspection:	In accordance with ISO 2859.	
Storage Instructions:	Store in a dry, cool place (<40°C) away from direct sunlight and fluorescent light.	
Shelf Life:	Five (5) years from date of manufacture.	
Bar Codes:	Each size carries a unique bar code on the inner dispenser box (EAN-13) and the outer carton (ITF-14). This allows automatic data entry into Stock Control Systems.	
Lot Number:	This appears on the sealed pouch, the inner dispenser box and the outer carton to provide traceability.	
Packaging:	One cuffed pair per inner wallet; One inner wallet per sealed peel pouch; 50 pouches per inner box; Four inner boxes per outer carton (200 pairs).	
Regulatory Classification:	Class 1 Sterile.	
Product Codes and Sizes:	ONBS201-S Small ONBS202-M Medium ONBS203-L Large ONBS204-XL X-Large	

ONBS/PDS7



GB04/61092