

Omega Techno™ Blue

NON-STERILE POWDERFREE NITRILE MEDICAL GLOVES



Master Code

Omega Techno™ Blue non-sterile powderfree ambidextrous nitrile medical gloves mark an important advance in non-latex thin glove technology.

These non-sterile powderfree ambidextrous nitrile medical gloves are 100% latex-free, eliminating allergy problems normally associated with natural latex proteins.

Our medical gloves have excellent chemical resistance and the soft formulation provides superior comfort and feel. They have an excellent grip and are available in a wide range of sizes to suit all users.



GB04/61092



ISO 9001
REGISTERED FIRM

NITRITEX LTD.

T: +44 (0) 1638 663338

E: info@nitritex.com

NITRITEX ASIA PTE LTD.

T: +65 6536 9123

E: nitritex@pacific.net.sg

NITRITEX PUERTO RICO Inc.

T: +1 787 309 3979

E: info@nitritex.com.pr

NITRITEX CANADA LTD.

T: +1 (905) 946 9539

E: info@nitritexcanada.com

NITRITEX MALAYSIA Sdn. Bhd.


T: +(603) 5569 3857

E: info@nitritex.com.my



Omega Techno™ Blue

NON-STERILE POWDERFREE NITRILE MEDICAL GLOVES

Colour:	Blue.	 Master Code
Length:	240mm (9.5").	
Material:	Nitrile.	
Type:	Ambidextrous.	
Surface:	Finger textured.	
Other Features:	Powderfree, beaded cuff.	
Protein Level:	Nitrile contains no natural latex proteins.	
Packing:	Boxed.	
AQL:	1.5 for holes (1000ml water test).	
Physical Properties:	Manufactured to exceed the requirements of the European Medical Gloves Standard EN455 (Parts 1, 2 & 3).	
Quality Systems:	Manufactured in a facility holding the ISO 9001:2000 registration.	
Sampling & Inspection:	In accordance with ISO 2859.	
Storage Instructions:	Store in a dry, cool place (<40°C) away from direct sunlight and fluorescent light.	
Shelf Life:	Five (5) years from date of manufacture.	
Bar Codes:	Each size carries a unique bar code on the inner dispenser box (EAN-13) and the outer carton (ITF-14). This allows automatic data entry into Stock Control Systems.	
Lot Number:	This appears on the inner dispenser box and the outer carton to provide traceability.	
Packaging:	100 pieces per inner dispenser box; 10 dispenser boxes per carton (1000 pieces).	
Re-order Codes and Sizes:	OTEBL2101-XS X-Small OTEBL2102-S Small OTEBL2103-M Medium OTEBL2104-L Large OTEBL2105-XL X-Large	

OTEBL/PD59



GB04/61092